

Deko-190 with an automatic 'Non-Touch' door



Deko-190 ward-disinfector is also available with an automatically opening and closing door.

An optical sensor activates the door mechanism and locking device. This eliminates the need to operate a door handle or foot pedal to open or close the door.

Simply move the article you are carrying (or your elbow or hand) across the front face of the optical sensor, which is fascia mounted. The door will automatically open/close for loading/unloading of items. An audible 'beep' will indicate the start of the movement.

If any obstruction occurs during the opening or closing operation of the door, it will stop the movement and reverse its operation.